

## REMEMBER!

**Always dispose of sharps at the points of use.  
Store containers safely with temporary closure  
across. Consider use of brackets.**

**Assemble correctly...snap lid on  
all around firmly and fill in label  
details**

**Place on bracket OR take to  
point of use in a tray.**

**Lock temporary  
closure when  
contents reach fill line**

**Ensure  
temporary  
closure is in the  
open position  
before use**

**Close temporary  
closure when not  
in use**

**Complete label on  
assembly and again  
prior to disposal**

**Remove container  
and dispose of  
according to your  
local policy**

